



**Hygienekontrolltechnik, Reinhold Pfeiffer**  
Chemnitzer Straße 33, D-09212 Limbach-Oberfrohna,  
Germany  
Tel./Fax +49 (0)3722-95370, E-mail: HK-Pfeiffer@web.de

## Chemical indicator for steam sterilization – PVP steam monitoring indicator

### Product information

<b>Field of application:</b>	Chemical steam sterilization indicator for batch control
<b>Benchmark standard:</b>	DIN EN ISO 11140-1, DIN EN ISO 15882
<b>Specifications:</b>	<p>Reaction behaviour: staged reaction, the colour changes from yellow via grey to black</p> <p>Reference: colour sample included as package insert for exposed saturated steam at 134°C with 3.5 min, Classification: visual Class 5, measure Class 6</p> <p>Indicator imprint: inverse image: STEAM and four squares</p> <p>Substrate: self-adhesive surface-protected paper, front side partially perforated and grooved</p> <p>Design: sheet with 25 indicator strips (75 mm x 6 mm), separately detachable, fits for the process challenge devices</p> <p>Packaging unit: 200 indicators on 8 sheets of 25 indicators each with a colour sample on the package insert</p> <p>Shelf life: non-sterilized indicators – 4 years with correct storage</p> <p>Storage: room temperature, protected from direct sunlight, relative humidity not more than 65 %</p> <p>Toxicological properties: no toxic substances are released</p>
<b>Instructions for use:</b>	<p>The indicator strips are inserted in the PCD and please note the product insert and the instructions of the Robert Koch Institute on. After sterilization, the indicators are compared with the colour sample and affixed to the documentation sheet by removing the anti-adhesive paper. The precise procedure is described in the package insert.</p> <p>The colour change is not reversible.</p>